

Quadrupole transport — 3D | Sets up a periodic FODO lattice and creates a beam that is matched to the lattice. | The beam is propagated one lattice period.

Outline

- Quadrupole transport 3D
- Quadrupole transport 2D-XY
- Plasma acceleration

Examples are in simbeamplasma18/warp_scripts from github repository:

• git clone https://github.com/uspas/simbeamplasma18.git.





FODO3D: tasks

- 1 Into FODO3D, open FODO3D.py
- 2 Execute file: "python -i FODO3D.py"
- ③ Open cgm file and explore:
 - a) "gist FOD03D.000.cgm &"
- ④ Compare the slices emittance diagnostics to the whole emittance diagnostic. Can you explain the differences (i.e. one is constant, the other oscillates)?
- ⑤ Read input script and try to understand every command
- 6 Change 'w3d.distrbtn = "semigaus" to 'w3d.distrbtn = "KV"; rerun & observe
- Thange 'w3d.distr_I = "gaussian" to 'w3d.distr_I = "neuffer"; rerun & observe
- (8) Insert "beam.x0 = beam.a0/2" on the line following "beam.a0 = ..."; rerun & observe
- 9 Check that you have the "ffmpeg" software installed: "which ffmpeg"
 - · If not, download and install ffmpeg
- ① Change "I_movieplot = False" to "I_movieplot = True" & rerun
 - If all goes well, after a few minutes, you should have a movie "movie.mp4"



Quadrupole transport – XY






File xy-quad-mag-mg.py:

nonrelativistic Warp xy slice simulation of a K⁺ ion beam with intense space-charge focused by a hard-edge magnetic quadrupole doublet focusing lattice.

Plasma acceleration



Scripts lpa_script.py, lpa_script_2d.py, pwfa_script.py - basic plasma acceleration runs:

 Generate plasma, laser or beam driver, and injected electron beam and follow self-consistent evolution.

xy-quad-mag-mg: tasks

- 1 Into XY-quad, open xy-quad-mag-mg.py
- ② Execute file: "python -i xy-quad-mag-mg.py"
- 3 Open cgm file and explore: "gist xy-quad-mag-mg.000.cgm &"
- 4 Read input script and try to understand every command.
- (5) Comment 'w3d.distrbtn = "SG"' and uncomment 'w3d.distrbtn = "KV"', rerun and compare to results using the KV vs SG distributions.
- ⑥ Change the initial emittance "emit = 10.e-6" to "emit = 10.e-7", rerun and observe effect on matching and emittance preservation.
- ① Change switch "I_automatch = False" to "I_automatch = True", rerun and observe difference with previous run.
- 8 Change n_grid=200 to 400, rerun and observe differences. Can you explain?
- With the simulation back at the python prompt, type 'dump()', then run for another 500 steps: "step(500)".
- 1 In another terminal, start python and type:
 - · from warp import *
 - restart('xy-quad-mag-mg001000.dump')
 - step(500)

Reopen "xy-quad-mag-mg.000.cgm" and compare to "xy-quad-mag-mg.001.cgm".

Plasma acceleration tasks

- ① In Plasma_acceleration, open LWFA_3D/lwfa_3d.py and PWFA_3D/pwfa_3d.py
- 2 Execute the files "python –i lwfa_3d.py" and "python –i pwfa_3d.py" separately.
- ③ It takes some time to run. While it runs, you may open periodically the cgm files ***_script.000.cgm and see the progress. In the meantime, also go through the input scripts and try to understand all the commands.
- 4 At the end of the run, a plot displays the energy of the accelerated beam versus z.
- (5) Explore data using the OpenPMD notebook viewer by running 'openPMD_notebook' in the run directory. While exploring data, run 2000 additional time steps.
- 6 In plasma_acceleration, open LWFA_2D/lwfa_2d.py and execute.
- ① Open cgm file and explore:
 - a) "gist lpa_basic_2d.000.cgm &"
- Read input script and try to understand every command.
- ① Run the script ptime displaying the history of the elapsed time and the time per step.
 Observe the spikes from the diagnostics.

